



INNOVATION IN SECURE DELIVERY

PREP-LOCK™

Tamper Evident Capping Station

Tamper-evident caps are an essential addition to any facilities' drug security protocols. With the Tamper Evident Capping Station, installing tamper-evident security could not be easier. This compact device automates cap installation and ensures proper application by providing perfect torque to every syringe, every time.

Productivity & Operational Efficiencies

Increase productivity and operational efficiencies with mechanized syringe capping.

Eliminates Pharmacist Fatigue

The Tamper Evident Capping Station eliminates pharmacist fatigue due to repeated capping.

Touch Free Capping

The touch-free capping process promotes aseptic technique.

Designed for ISO 5 & GMP Requirements

The easy-to-clean design is thoughtfully crafted for use in ISO 5 environments and allows GMP requirements to be conducted efficiently and safely.



MADE IN THE USA

IMIWEB.COM | 1.800.344.2554



Prep-Lock™ | Tamper Evident Capping Station

| Order Number | Description | Qty |
|--------------|-----------------|--------|
| TECS01 | Capping Station | 1 Unit |

| Technical Factors | |
|--------------------|--|
| Latex: | •This product is not made with natural rubber latex. |
| Product Material: | •Made with anodized aluminum, stainless steel & delrin. |
| Conflict Material: | •IMI products contain no conflict minerals. |
| Cleaning: | •Compatible with 70% IPA, Decon-Quat®/ Decon-PHene® and Peridox-RTU® cleaning agents. |
| Certification: | ETL Certified Conforms to UL Standard 61010-1 Certified to CSA Standard C22.2#61010-1-12 |
| Warranty: | 1 Year limited warranty |

To place an order, please call or visit

IMIWEB.COM | 1.800.344.2554

