



Prep-Lock™

Tamper Evident Caps

For ENFit® Syringes

FEATURING RFID

Prevent Diversion Events

Actively deter diversion by providing a physical barrier and powerful analytic capabilities to your drugs at product level.

Traceability & Visibility

DoseID RFID takes serialization from the unit of sale to the unit of use. It uniquely serializes every dose to track each action taken upon it during its inventory life cycle.

Minimal staff investment or procedural changes.

Pre-tagged Tamper Evident Caps can reduce the cost and time of RFID implementation.

Simplify Labor Intensive Processes

Improve ordering and supply chain management with unit-of-use medication serialization.



MADE IN THE USA

IMIWEB.COM | 1.800.344.2554



IMI is a proud member of DoseID, a consortium of industry innovators established to unify, develop and administer comprehensive RFID standards for the pharmaceutical Healthcare industry.



IMI-501-AA-71 R1



Prep-Lock™ | Tamper Evident Caps For ENFit Syringes FEATURING RFID

Order Number	Description	Qty
57-220-100 RFID	For use with ENFit® Syringes Purple	Box of 100

Technical Factors	
Sterility:	•Packaged non sterile.
DEHP:	•This product does not contain DEHP.
BPA:	•This product does not contain BPA.
Latex:	•This product is not made with natural rubber latex.
Product Material:	•Made with polypropylene & polystyrene.
DoseID RFID:	<ul style="list-style-type: none"> •Products carrying DoseID certification ensure the quality, performance, and interoperability of RFID tagged products with DoseID certified systems. •DoseID certified RFID tags are interoperable with GTIN, SGTIN, DSCSA, block chain & other standards. •Dose ID tags are compliant with existing RAIN and GSI standards.

ENFit® is a GEDSA trademark.

To place an order, please call or visit
IMIWEB.COM | 1.800.344.2554

MADE IN THE
USA

