



Prep-Lock™

# Tamper Evident Caps For Oral Syringes

**FEATURING RFID**

## Prevent Diversion Events

Actively deter diversion by providing a physical barrier and powerful analytic capabilities to your drugs at product level.

## Traceability & Visibility

DoseID RFID takes serialization from the unit of sale to the unit of use. It uniquely serializes every dose to track each action taken upon it during its inventory life cycle.

## Minimal staff investment or procedural changes.

Pre-tagged Tamper Evident Caps can reduce the cost and time of RFID implementation.

## Simplify Labor Intensive Processes

Improve ordering and supply chain management with unit-of-use medication serialization.



IMI is a proud member of DoseID, a consortium of industry innovators established to unify, develop and administer comprehensive RFID standards for the pharmaceutical Healthcare industry.



**MADE IN THE USA**

IMIWEB.COM | 1.800.344.2554



# Prep-Lock™ | Tamper Evident Caps For Oral Syringes FEATURING RFID

Each Oral Syringe Cap is designed to the unique syringe manufacture specifications.

Order Number	Description	Qty
57-70-100 RFID	Compatible with Baxter Exacta-Med Oral Syringes Red Top, Non-Sterile	Box of 100
57-73-100 RFID	Compatible with BD Univia™ Oral Syringes Yellow Top, Non-Sterile	Box of 100

Technical Factors	
Sterility:	•Packaged non sterile.
DEHP:	•This product does not contain DEHP.
BPA:	•This product does not contain BPA.
Latex:	•This product is not made with natural rubber latex.
Product Material:	•Made with polypropylene & polystyrene.
DoseID RFID:	•Products carrying DoseID certification ensure the quality, performance, and interoperability of RFID tagged products with DoseID certified systems. •DoseID certified RFID tags are interoperable with GTIN, SGTIN, DSCSA, block chain & other standards. •Dose ID tags are compliant with existing RAIN and GSI standards.

To place an order, please call or visit

IMIWEB.COM | 1.800.344.2554

MADE IN THE  
**USA**



IMI-501-AA-72 R1